

Nippon Protein Co., Ltd. Product Specification	Established Date	Apr. 2, 2012	Spec. No.	URA-10055N-1	Page	1/4
	Effective Date	Mar. 12, 2021				

**N-Acetyl-L-Cysteine(EX) \*1**  
: For Manufacturing, Processing or Repacking

**N-Ac-L-Cys(EX)**

C<sub>5</sub>H<sub>9</sub>NO<sub>3</sub>S:163.19

N-Acetyl-L-Cysteine(EX), when dried, contains not less than 98.5 percent and not more than 101.0 percent of N-Acetyl-L-Cysteine (C<sub>5</sub>H<sub>9</sub>NO<sub>3</sub>S).

**Description:** White crystals or crystalline powder; strongly acid taste. Freely soluble in water, soluble in ethanol (96).

**Identification:** Compare the infrared absorption spectrum <NP TEST 35> of the sample with that of the standard by potassium bromide disc method.

**Specifications:**

Item	Limit	Test Method
<b>Specific rotation</b> [ $\alpha$ ] <sup>20</sup> <sub>D</sub>	+21.3 to +27.0°	NP TEST 1 [dried sample, C=5, phosphate buffer solution*2]*3
<b>State of solution</b> (Transmittance)	clear and colorless not less than 98.0%	NP TEST 2 [0.5g in 10mL of H <sub>2</sub> O, spectrophotometer, 430nm, 10mm cell thickness]
<b>Heavy metals(Pb)</b>	not more than 10ppm	NP TEST 7 [1.0g, (3), ref: 1.0mL of Pb Std. (0.01mg/mL)]
<b>Zinc(Zn)</b>	not more than 10.0ppm	EP
<b>Related substances</b>	not more than 0.2% of L-cysteine not more than 0.3% of N, N'-diacetyl-L-cystine not more than 0.15% of N, S-diacetyl-L-cysteine not more than 0.10% of other impurity (each) not more than 0.5% of total impurities	EP*4
<b>Loss on drying</b>	not more than 1.00%	NP TEST 11 [1g, in vacuum, at 70°C for 4 hours]

Nippon Protein Co., Ltd. Product Specification	Established Date	Apr. 2, 2012	Spec. No.	URA-10055N-1	Page	2/4
	Effective Date	Mar. 12, 2021				

Item	Limit	Test Method
<b>Residue on ignition</b> (sulfated)	not more than 0.20%	NP TEST 13 [1g, at 550°C to 650°C for 3 hours]
<b>Assay</b>	98.5 to 101.0%	EP
<b>pH</b>	2.0 to 2.8	NP TEST 33 [1.0g in 100mL of H <sub>2</sub> O]

<b>Endotoxin</b>	less than 24.0EU/g	NP TEST 34 [C = 0.25, kinetic-turbidimetric technique]
------------------	--------------------	--

\* 1 : USP, EP

\* 2 : In a 50mL volumetric flask mix 2.5g of the test sample, accurately weighed, with 2mL of disodium dihydrogen ethylenediaminetetraacetate solution (1→100), add 15mL of sodium hydroxide solution (1→25), and mix to dissolve. Dilute to volume with pH 7.0 buffer prepared by mixing 29.5mL of 1mol/L sodium hydroxide, 50mL of 1mol/L potassium dihydrogen phosphate, and sufficient water to make 100mL.

\* 3 : Temperature coefficient of  $[\alpha]_D^{25}$  :  $-0.06^\circ$

\* 4 : 1) Analytical condition

Detector: UV detector (220nm)

Column: ODS (4.0mm  $\phi$  × 0.25m)

Temperature: Constant temperature around 25°C

Mobile phase: Acetonitrile, water previously adjusted to pH3.0 with phosphoric acid (3:97 V/V)

Flow rate: 1.0mL/min

Injection volume: 20 $\mu$ L

Analytical time: Three times over the retention time for acetylcysteine

2) Solution preparation

i ) Test solution

Take 800mg of sample and add 0.01mol/L hydrochloric acid solution to make 100mL.

ii ) Reference solution(a)

Dilute 5.0mL of the test solution to 50.0mL with 0.01mol/L hydrochloric acid solution. Dilute 1.0mL of this solution to 100.0mL with 0.01mol/L hydrochloric acid solution.

iii) Reference solution (b)

Dissolve 80mg of L-cystine in 2mL of 1mol/L hydrochloric acid solution and dilute to 200mL with water.

Nippon Protein Co., Ltd. Product Specification	Established Date	Apr. 2, 2012	Spec. No.	URA-10055N-1	Page	3/4
	Effective Date	Mar. 12, 2021				

iv) Reference solution (c)

Dissolve 100mg of L-cysteine in 100mL of 0.01mol/L hydrochloric acid solution (L-cysteine dilute solution).

Dissolve 5mg of N,N'-diacetyl-L-cystine and 2.5mg of N,S-diacetyl-L-cysteine in 0.01mol/L hydrochloric acid solution, mix with 3mL of L-cysteine dilute solution and 4mL of Reference solution (b) and dilute to 20mL with 0.01mol/L hydrochloric acid solution. Dilute 5mL of this solution to 50 mL with the test solution.

v) Reference solution (d)

Dissolve 4mg of sodium 2-methyl-2-thiazoline-4-carboxylate in 0.01mol/L hydrochloric acid solution and dilute to 100mL with 0.01mol/L hydrochloric acid solution.

3) Calculation

Impurity contents (percentage) are calculated according to following equations.

i ) In case of L-cysteine, N, N'-diacetyl-L-cystine, N, S'-diacetyl-L-cysteine

$$\text{Each impurity (\%)} = (A1 / A2) \times (m2 / m1) \times 100 \times \text{correction factor}$$

ii ) In case of unspecified impurities

$$\text{Each impurity (\%)} = (A3 / A2) \times (m2 / m1) \times 100$$

Where,

A<sub>1</sub>: peak area of individual impurity (L-cysteine, N,N'-diacetyl-L-cystine and N,S'-diacetyl-L-cysteine) in the chromatogram obtained with test solution.

A<sub>2</sub>: peak area of acetylcysteine in the chromatogram obtained with reference solution (a).

A<sub>3</sub>: peak area of unspecified impurity in the chromatogram obtained with test solution.

m<sub>1</sub>: mass of the acetylcysteine in test solution (mg/mL)

m<sub>2</sub>: mass of acetylcysteine in reference solution (a) (mg/mL)

correction factor: L-cysteine = 3.4

N,N'-diacetyl-L-cystine = 0.7

N,S'-diacetyl-L-cysteine = 0.3

Nippon Protein Co., Ltd. Product Specification	Established Date	Apr. 2, 2012	Spec. No.	URA-10055N-1	Page	4/4
	Effective Date	Mar. 12, 2021				

(Remark)

Following peaks are neglected;

- i ) Any peak originated from solvent used for mobile phase.
- ii ) Peak due to 2-methyl-2-thiazoline-4-carboxylic acid, which is formed due to degradation of acetylcysteine in acidic solution.

End of document