L-Arginine Monohydrochloride¹

C₆H₁₄N₄O₂·HCl: 210.66

L-Arginine Monohydrochloride, when dried, contains not less than 98.5 percent and not more than 101.0 percent of L-Arginine Monohydrochloride ($C_6H_{14}N_4O_2$ ·HCl).

Description

White crystals or crystalline powder; slightly characteristic taste. Freely soluble in water and in formic acid, very slightly soluble in ethanol (95).

Identification

Compare the infrared absorption spectrum of the sample with that of the standard by potassium bromide disc method.

Specifications

Item	Limit	Test
Specific rotation $\left[\alpha\right]_{D}^{20}$	+22.1 to +22.9°	AJI TEST 1
		[Dried sample, C=8, 6mol/L HCl] ²
State of solution (Transmittance)	Clear and colorless	AJI TEST 2
	Not less than 98.0%	[1.0g in 10mL of H ₂ O, spectrophotometer, 430nm, 10mm cell
		thickness]
Chloride (Cl)	16.58 to 17.00%	AJI TEST 3
		[Dried sample, 1.0g ³ , B]
Ammonium (NH ₄)	Not more than 0.02%	AJI TEST 4
		[D-1]
Sulfate (SO ₄)	Not more than 0.020%	AJI TEST 5
		[0.85g, (1), ref: 0.35mL of 0.005mol/L H ₂ SO ₄]
Iron (Fe)	Not more than 10ppm	AJI TEST 6
		[0.75g, B-1, ref: 0.75mL of Iron Std. (0.01mg/mL)]
Heavy metals (Pb)	Not more than 10ppm	AJI TEST 7
		[1.0g, (1), ref: 1.0mL of Pb Std. (0.01mg/mL)]
Arsenic (As ₂ O ₃)	Not more than 1ppm	AJI TEST 8
		[2.0g, (1), ref: 2.0mL of As ₂ O ₃ Std.]
Related substances	1)Conforms ⁴	AJI TEST 9
		[Test sample: 50µg, S-6-a, control; L-Arg·HCl 0.25µg]
	2)Any unspecified impurity	AJI TEST 26 ⁵
	Not more than 0.20%	
	Total impurities	
	Not more than 0.50%	
Loss on drying	Not more than 0.20%	AJI TEST 11
		[1g, at 105°C for 3 hours]
Residue on ignition	Not more than 0.10%	AJI TEST 13
(Sulfated)		[1g, at 550°C to 650°C for 3 hours]

Amino Acids Specifications / Monographs		
L-Arginine Monohydrochloride		2/2
Issued Date: Jun. 28, 2019		

L-Arginine Monohydrochloride¹

Specifications (cont'd)

Item	Limit	Test
Assay	98.5 to 101.0%	AJI TEST 14
		[Dried sample, 100mg, (3), 2mL of formic acid, 0.1mol/L
		HClO ₄ 1mL=10.53mg C ₆ H ₁₄ N ₄ O ₂ ·HCl]
рН	4.7 to 6.2	AJI TEST 33
		[1.0g in 10mL of H ₂ O]

The test for Endotoxin when the material will be used for manufacturing parenteral products is as follows:

Item	Limit	Test
Endotoxin	Less than 6.0EU/g	AJI TEST 34
		[C=1, kinetic-turbidimetric technique]

¹ This product, in terms of actual quality, conforms to USP, EP, and JP.

² Temperature coefficient of $[\alpha]_{D}^{t}$: -0.03°

⁵ Disregard limit: 0.05%

End of document

³ Weigh accurately 1.0g of the sample and add water to make exactly 100mL.

⁴ Any secondary spot in the chromatogram obtained from the Test Solution is less intense than the principal spot in the chromatogram obtained from the Standard Solution: the number of those spots is not more than four and not more than 2.0% of the total impurities found.